



Rapidfix

Dual Lead Thoracolumbar Screws



These screws are used with
Tower Fix Scoliosis
Deformity Instrument Set
and Plexus & Silver
Posterior Spinal Instrument
Set

(Set Code: N-SS-002, Set Code: N-SS-001)



Rapidfix Dual Lead Polyaxial Screw (5 mm increments)

Product Code	Diameter(mm)	XX/L(mm)
NSPAK45XX	4.5	20-45
NSPAK50XX	5.0	25-50
NSPAK55XX	5.5	20-55
NSPAK60XX	6.0	30-55
NSPAK65XX	6.5	25-60
NSPAK70XX	7.0	30-55
NSPAK75XX	7.5	25-100
NSPAK80XX	8.0	25-100



Rapidfix Dual Lead Reduction Polyaxial Screw (5 mm increments)

Product Code	Diameter(mm)	XX/L(mm)
NSPAKR45XX	4.5	25-45
NSPAKR50XX	5.0	25-50
NSPAKR55XX	5.5	20-55
NSPAKR60XX	6.0	30-55
NSPAKR65XX	6.5	25-60
NSPAKR70XX	7.0	30-55
NSPAKR75XX	7.5	25-100
NSPAKR80XX	8.0	25-100



L: Length(mm)

Features;

- Rapidfix Dual Lead screw design displays a superior pullout compared to competitor design fast and controlled Insertion
- Combination of two threads having the same form of a double threaded.
- Excellent tactile feedback on bone purchase of screw thanks to combination of double lead

Rapidfix Dual Lead Cannulated Fenestrated Polyaxial Screw (5 mm increments)
Product Code D(mm) XX/L(mm)

NSPA KH55XX	5.5	30-55
NSPA KH60XX	6.0	30-55
NSPA KH65XX	6.5	30-55
NSPA KH70XX	7.0	30-55
NSPA KH75XX	7.5	30-55
NSPA KH80XX	8.0	30-55



Rapidfix Dual Lead Monoaxial Screw (5 mm increments)
Product Code D(mm) XX/L(mm)

N SMAK40XX	4.0	20-45
N SMAK45XX	4.5	20-45
N SMAK50XX	5.0	25-50
N SMAK55XX	5.5	20-55
N SMAK60XX	6.0	30-55
N SMAK65XX	6.5	25-60
N SMAK70XX	7.0	30-55
N SMAK75XX	7.5	25-100
N SMAK80XX	8.0	25-100



Rapidfix Dual Lead Reduction Monoaxial Screw (5 mm increments)
Product Code D(mm) XX/L(mm)

N SMAKR45XX	4.5	20-45
N SMAKR50XX	5.0	25-50
N SMAKR55XX	5.5	20-55
N SMAKR60XX	6.0	30-55
N SMAKR65XX	6.5	25-60
N SMAKR70XX	7.0	30-55
N SMAKR75XX	7.5	25-100
N SMAKR80XX	8.0	25-100



Plexus Rod
Product Code D(mm) L(mm) Piece

NSR15504	5.5	40	
NSR15505	5.5	50	
NSR15506	5.5	60	
NSR15507	5.5	70	
NSR15508	5.5	80	
NSR15509	5.5	90	
NSR15510	5.5	100	
NSR15511	5.5	110	
NSR15512	5.5	120	
NSR15513	5.5	130	
NSR15514	5.5	140	
NSR15540	5.5	400	
NSR15550	5.5	500	
NSR15560	5.5	600	
NSR1604	6.0	40	
NSR1605	6.0	50	
NSR1606	6.0	60	
NSR1607	6.0	70	
NSR1608	6.0	80	
NSR1609	6.0	90	
NSR1610	6.0	100	
NSR1611	6.0	110	
NSR1612	6.0	120	
NSR1613	6.0	130	
NSR1614	6.0	140	
NSR1615	6.0	150	
NSR1616	6.0	160	
NSR1617	6.0	170	
NSR1618	6.0	180	
NSR1619	6.0	190	
NSR1620	6.0	200	
NSR1621	6.0	210	
NSR1622	6.0	220	




Rod
Product Code D(mm) L(mm) Piece

NSR1623	6.0	230	
NSR1624	6.0	240	
NSR1625	6.0	250	
NSR1626	6.0	260	
NSR1627	6.0	270	
NSR1628	6.0	280	
NSR1629	6.0	290	
NSR1630	6.0	300	
NSR1631	6.0	310	
NSR1632	6.0	320	
NSR1640	6.0	400	
NSR1648	6.0	480	
NSR1650	6.0	500	
NSR1660	6.0	600	



Cr-Co Rod
Product Code D(mm) L(mm) Piece

N SCRR15504	5.5	40	
N SCRR15505	5.5	50	
N SCRR15506	5.5	60	
N SCRR15507	5.5	70	
N SCRR15508	5.5	80	
N SCRR15509	5.5	90	
N SCRR15510	5.5	100	
N SCRR15511	5.5	110	
N SCRR15512	5.5	120	
N SCRR15513	5.5	130	
N SCRR15514	5.5	140	
N SCRR15540	5.5	400	
N SCRR15550	5.5	500	
N SCRR15560	5.5	600	
N SCRR1604	6.0	40	
N SCRR1605	6.0	50	



D: Diameter(mm), L: Length(mm)

Cannulated and Fenestrated screws are used with a 3.5 hex screwdriver.



This implant system used with **Tower Fix Scoliosis Deformity Instrument Set** and **Plexus & Silver Posterior Spinal Instrument Set**
 (Set Code: N-SS-002, Set Code: N-SS-001)

Rapidfix
 Dual Lead Thracolumbar
 Screws